MSDS

Material Safety Data Sheet

Date prepared: Sep 5, 2012 Date revised: N/A (first edition)

Section 1 - General Information

Chemical Name: Methacrolein 2,4-dinitrophenylhydrazone

CAS Reg. Number: 5077-73-6
Catalog Number: 1ST1027
Manufacturer's Name: A Chemtek, Inc.

100 Barber Avenue

Worcester, MA 01606, USA

Telephone Number: 508-471-4121 **Fax Number:** 508-845-9201

For R & D use only

Section 2 - Hazardous Ingredients/Identity Information

Classification according to Regulation (EC) No 1272/2008 [EU-GHS/CLP]

Skin irritation (Category 2) Eye irritation (Category 2)

Specific target organ toxicity - single exposure (Category 3)

Classification according to EU Directives 67/548/EEC or 1999/45/EC

Irritating to eyes, respiratory system and skin.

Label elements

Labelling according Regulation (EC) No 1272/2008 [CLP]

!>

Pictogram
Signal word
Warning

Potential Health Effects

Inhalation: May be harmful if inhaled. May cause respiratory tract irritation. Skin: May be harmful if absorbed through skin. May cause skin irritation.

Eyes: May cause eye irritation.
Ingestion: May be harmful if swallowed..

Section 3 - Composition/Information on Ingredients

Formula: $C_{10}H_{10}N_4O_4$ Molecular Weight: 250.22g/molCAS-No. 5077-73-6

Section 4 - First Aid Measures

If inhaled

If breathed in, move person into fresh air. If not breathing give artificial respiration.

In case of skin contact

Wash off with soap and plenty of water.

In case of eye contact

Flush eyes with water as a precaution.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

Section 5 - Fire-Fighting Measures

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for fire-fighters

Wear self-contained breathing apparatus for firefighting if necessary.

Section 6 - Accidental Release Measures

Personal precautions Avoid dust formation

Environmental precautions Do not let product enter drains.

Methods for cleaning up Sweep up and shovel. Keep in suitable, closed containers for disposal.

Section 7 - Handling and Storage

MSDS

Material Safety Data Sheet

Date prepared: Sep 5, 2012 Date revised: N/A (first edition)

Handling Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for

preventive fire protection. Avoid contact with skin and eyes. Avoid formation of dust and

aerosols.

Storage Keep container tightly closed in a dry and well-ventilated place.

Section 8 - Exposure Controls/Personal Protection

Contains no substances with occupational exposure limit values

Personal protective equipment

Respiratory protection

For nuisance exposures use type P95 (US) or type P1 (EU EN 143) particle respirator. For higher level protection use type OV/AG/P99 (US) or type ABEK-P2 (EU EN 143) respirator cartridges. Use Aldrich - 429678 Page 4 of 6 respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection For prolonged or repeated contact use protective gloves

Eye protection Safety glasses

Hygiene measures General industrial hygiene practice

Section 9 - Physical and Chemical Properties

Appearance:

Form Solid Color Orange-red

Safety data

pH: not available Melting point: 204-206°C - lit. Boiling point: not available Flash point: not available Ignition temperature: not available Lower explosion limit: not available Upper explosion limit: not available Water solubility: not available

Section 10 - Stability and Reactivity

Storage stability Stable under recommended storage conditions

Materials to avoid Strong oxidizing agents

Hazardous decomposition products Other decomposition products - no data available

Section 11 - Toxicological Information

Acute toxicity:not availableIrritation and corrosion:not availableSensitization:not availableRespiratory or skin sensitization:not available

Chronic exposure

No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

Potential health effects

Inhalation: May be harmful if inhaled. May cause respiratory tract irritation. Skin: May be harmful if absorbed through skin. May cause skin irritation.

Eyes: May cause eye irritation.
Ingestion: May be harmful if swallowed.

Section 12 - Ecological Information

Elimination information (persistence and degradability)

not available

Ecotoxicity effects: not available **Further information on ecology:** not available

Section 13 - Disposal Considerations

MSDS

Material Safety Data Sheet

Date prepared: Sep 5, 2012 Date revised: N/A (first edition)

Product Observe all federal, state, and local environmental regulations.

Contaminated packaging Dispose of as unused product.

Section 14- Transport Information

DOT (US):

IMDG:

Not dangerous goods

Not dangerous goods

Not dangerous goods

Section 15 - Regulatory Information

This safety datasheet complies with the requirements of Regulation (EC) No. 1907/2006.

Section 16 - Other Information

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. A Chemtek shall not be held liable for any damage resulting from handling or from contact with the above product.