

MSDS

Material Safety Data Sheet

Date prepared: August, 29th, 2012

Date revised: N/A (first edition)

Section 1 - General Information

Chemical Name: Dicyclohexyl phthalate
Synonyms: NA
CAS Reg. Number: 84-61-7
Catalog Number: 1ST1115
Manufacturer's Name: A Chemtek, Inc.
100 Barber Avenue
Worcester, MA 01606, USA
Telephone Number: 508-471-4121
Fax Number: 508-845-9201

For R & D use only

Section 2 - Hazardous Ingredients/Identity Information

Emergency Overview

OSHA Hazards

Irritant

GHS Classification

Skin irritation (Category 2)

Eye irritation (Category 2A)

Specific target organ toxicity - single exposure (Category 3)

GHS Label elements, including precautionary statements



Pictogram

Signal word

Warning

Hazard statement(s)

H315 Causes skin irritation.
H319 Causes serious eye irritation.
H335 May cause respiratory irritation.

Precautionary statement(s)

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray
P264 Wash skin thoroughly after handling.
P271 Avoid release to the environment.
P280 Wear protective gloves/ eye protection/ face protection.
P302 + P352 IF ON SKIN: Wash with plenty of soap and water
P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P302 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P312 Call a POISON CENTER or doctor/ physician if you feel unwell.
P321 Specific treatment (see supplemental first aid instructions on this label).
P405 Store locked up.
P501 Dispose of contents/ container to an approved waste disposal plant.

HMIS Classification

Health hazard: 2
Chronic Health Hazard: *
Flammability: 0
Physical hazards: 0

NFPA Rating

Health hazard: 2
Fire: 0
Reactivity Hazard: 0

Section 3 - Composition/Information on Ingredients

Formula: $C_{20}H_{26}O_4$
Molecular Weight: 330.42 g/mol

MSDS

Material Safety Data Sheet

Date prepared: August, 29th, 2012

Date revised: N/A (first edition)

CAS-No. 84-61-7

Section 4 - First Aid Measures

GENERAL ADVICE: Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

EYE CONTACT: Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

SKIN CONTACT: Wash off with soap and plenty of water. Consult a physician.

INHALATION: If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

INGESTION: Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

Section 5 - Fire-Fighting Measures

Flammable properties

Flash point not available
Ignition temperature not available

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for fire-fighters

Wear self-contained breathing apparatus for firefighting if necessary.

Section 6 - Accidental Release Measures

Personal precautions Use personal protective equipment. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas.

Environmental precautions Prevent further leakage or spillage if safe to do so. Do not let product enter drains.

Methods for cleaning up Sweep up and shovel. Keep in suitable, closed containers for disposal.

Section 7 - Handling and Storage

Handling Avoid inhalation of vapour or mist.

Storage Keep container tightly closed in a dry and well-ventilated place.

Section 8 - Exposure Controls/Personal Protection

Components with workplace control parameters

Personal protective equipment

Respiratory protection Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Eye protection Safety glasses with side-shields conforming to EN166 Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin protection Impervious clothing, The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Hygiene measures Handle in accordance with good industrial hygiene and safety practice. Wash hands before breaks and at the end of workday.

MSDS

Material Safety Data Sheet

Date prepared: August, 29th, 2012

Date revised: N/A (first edition)

Section 9 - Physical and Chemical Properties

Appearance: Form Soild
Color not available

Safety data

pH: not available
Melting point: Melting point/range: 63 - 67 °C (145 - 153 °F) - lit.
Boiling point: not available
Flash point: not available
Ignition temperature: not available
Lower explosion limit: not available
Upper explosion limit: not available
Water solubility: not available

Section 10 - Stability and Reactivity

Storage stability Stable under recommended storage conditions.

Materials to avoid Strong oxidizing agents.

Hazardous decomposition products Hazardous decomposition products formed under fire conditions.
Carbon oxides.

Section 11 - Toxicological Information

Acute toxicity:

Oral LD50

not available

Inhalation LC50

not available

Dermal LD50

not available

Other information on acute toxicity

not available

Irritation and corrosion: not available

Sensitization: not available

Respiratory or skin sensitization: not available

Chronic exposure No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC, ACGIH, NTP, or OSHA.

Section 12 - Ecological Information

Elimination information (persistence and degradability)

not available

Ecotoxicity effects: not available

Further information on ecology: not available

Section 13 - Disposal Considerations

Product Offer surplus and non-recyclable solutions to a licensed disposal company.
Contact a licensed professional waste disposal service to dispose of this material.

Contaminated packaging: Dispose of as unused product.

Section 14- Transport Information

DOT (US):

Not dangerous goods

IMDG:

Not dangerous goods

IATA:

MSDS

Material Safety Data Sheet

Date prepared: August, 29th, 2012

Date revised: N/A (first edition)

Not dangerous goods

Section 15 - Regulatory Information

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

Acute Health Hazard

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

Dicyclohexyl phthalate	CAS-No. 84-61-7	Revision Date	2009-07-17
------------------------	-----------------	---------------	------------

New Jersey Right To Know Components

Dicyclohexyl phthalate	CAS-No. 84-61-7	Revision Date	2009-07-17
------------------------	-----------------	---------------	------------

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Section 16 - Other Information

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. A Chemtek shall not be held liable for any damage resulting from handling or from contact with the above product.