Material Safety Data Sheet

Date prepared: Oct 25, 2015 Date revised: N/A (first edition)

Section 1 - General Information

Chemical Name:FlumethrinCAS Reg. Number:69770-45-2Catalog Number:1ST20238Manufacturer's Name:A Chemtek, Inc.

100 Barber Avenue

Worcester, MA 01606, USA

Telephone Number: 508-471-4121 **Fax Number:** 508-845-9201

For R & D use only

Section 2 - Hazardous Ingredients/Identity Information

Classification of the substance or mixture

Acute toxicity, Oral (Category 3), H301 Acute toxicity, Inhalation (Category 3), H331 Acute toxicity, Dermal (Category 3), H311 Skin irritation (Category 2), H315 Eye irritation (Category 2A), H319 Skin sensitisation (Category 1), H317



Pictogram
Signal word
Danger

Hazard statement(s)

H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.
H319 Causes serious eye irritation.

Precautionary statement(s)

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P264 Wash skin thoroughly after handling.

P270 Do not eat, drink or smoke when using this product.
P271 Use only outdoors or in a well-ventilated area.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves/ protective clothing.

P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER or doctor/

physician.

P302 + P350 IF ON SKIN: Gently wash with plenty of soap and water.

P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position

comfortable for breathing.

P310 Immediately call a POISON CENTER or doctor/ physician.

P322 Specific measures (see supplemental first aid instructions on this label).

P330 Rinse mouth.

P361 Remove/Take off immediately all contaminated clothing.

P363 Wash contaminated clothing before reuse.

P391 Collect spillage.

P403 + P233 Store in a well-ventilated place. Keep container tightly closed.

P405 Store locked up.

P501 Dispose of contents/ container to an approved waste disposal plant.

Section 3 - Composition/Information on Ingredients

Material Safety Data Sheet

Date prepared: Oct 25, 2015 Date revised: N/A (first edition)

Formula: C28H22Cl2FNO3 Molecular Weight: 510.38g/mol CAS-No. 69770-45-2

Section 4 - First Aid Measures

If inhaled

If breathed in, move person into fresh air. If not breathing give artificial respiration.

In case of skin contact

Wash off with soap and plenty of water.

In case of eye contact

Flush eyes with water as a precaution.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

Section 5 - Fire-Fighting Measures

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for fire-fighters

Wear self-contained breathing apparatus for firefighting if necessary.

Section 6 - Accidental Release Measures

Personal precautions Avoid dust formation. Avoid breathing vapours, mist or gas.

Environmental precautions Do not let product enter drains.

Methods for cleaning up Sweep up and shovel. Keep in suitable, closed containers for disposal.

Section 7 - Handling and Storage

Handling Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for

preventive fire protection.

Storage Keep container tightly closed in a dry and well-ventilated place. Containers which are opened

must be carefully resealed and kept upright to prevent leakage.

Recommended storage temperature: 2 - 8 °C

Section 8 - Exposure Controls/Personal Protection

Contains no substances with occupational exposure limit values

Personal protective equipment

Respiratory protection

Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multipurpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection For prolonged or repeated contact use protective gloves

Eye protection Safety glasses

Hygiene measures General industrial hygiene practice

Section 9 - Physical and Chemical Properties

Appearance:

Form liquid

Color not available

Safety data

pH: not available
Melting point: not available
Boiling point: not available

Material Safety Data Sheet

Date prepared: Oct 25, 2015 Date revised: N/A (first edition)

Flash point: not available
Lower explosion limit: not available
Upper explosion limit: not available
Water solubility: not available

Section 10 - Stability and Reactivity

Storage stability Stable under recommended storage conditions

Materials to avoid Strong oxidizing, agents Strong bases, Strong oxidizing agents
Hazardous decomposition products
Hazardous decomposition products formed under fire conditions.

- Carbon oxides, nitrogen oxides (NOx), Hydrogen chloride gas, Hydrogen fluoride.

Section 11 - Toxicological Information

Acute toxicity:not availableIrritation and corrosion:not availableSensitization:not availableRespiratory or skin sensitization:not available

Chronic exposure

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

Potential health effects

Inhalation: not available
Skin: not available
Eyes: not available
Ingestion: not available

Section 12 - Ecological Information

Elimination information (persistence and degradability) not available

Ecotoxicity effects: not available **Further information on ecology:** not available

Section 13 - Disposal Considerations

Product Offer surplus and non-recyclable solutions to a licensed disposal company.

Contaminated packaging Dispose of as unused product.

Section 14- Transport Information

DOT (US):

UN number: 2810 Class: 6.1 Packing group: II

Proper shipping name: Toxic, liquids, organic, n.o.s. (α-Cyano-4-fluoro-3-phenoxybenzyl 3-[2-chloro-2-(4-

chlorophenyl)vinyl]-2,2-dimethylcyclopropanecarboxylate)

Marine pollutant: No

Poison Inhalation Hazard: No

IMDG:

UN number: 2810 Class: 6.1 Packing group: II EMS-No: F-A, S-A

Proper shipping name: TOXIC LIQUID, ORGANIC, N.O.S. (α-Cyano-4-fluoro-3-phenoxybenzyl 3-[2-chloro-2-(4-

chlorophenyl)vinyl]-2,2-dimethylcyclopropanecarboxylate)

Material Safety Data Sheet

Date prepared: Oct 25, 2015 Date revised: N/A (first edition)

Marine pollutant: No

IATA:

UN number: 2810 Class: 6.1 Packing group: II

Proper shipping name: Toxic liquid, organic, n.o.s. (α-Cyano-4-fluoro-3-phenoxybenzyl 3-[2-chloro-2-(4-

chlorophenyl)vinyl]-2,2-dimethylcyclopropanecarboxylate)

Section 15 - Regulatory Information

SARA 302 Components

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

α-Cyano-4-fluoro-3-phenoxybenzyl 3-[2-chloro-2-(4- chlorophenyl)vinyl]-2,2-dimethylcyclopropanecarboxylate CAS-No 69770-45-2 Revision Date

New Jersey Right To Know Components

α-Cyano-4-fluoro-3-phenoxybenzyl 3-[2-chloro-2-(4- chlorophenyl)vinyl]-2,2-dimethylcyclopropanecarboxylate CAS-No 69770-45-2 Revision Date

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Section 16 - Other Information

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. A Chemtek shall not be held liable for any damage resulting from handling or from contact with the above product.