MSDS

Date prepared: July 25, 2014

Material Safety Data Sheet

Date revised: N/A (first edition)

Section 1 - General Information

Chemical Name: Synonyms: CAS Reg. Number: Catalog Number: Manufacturer's Name: Melengestrol acetatesolution , 100ppm NA 2919-66-6 1ST2209-100A A Chemtek, Inc. 100 Barber Avenue Worcester, MA 01606, USA 508-471-4121 508-845-9201

Telephone Number: Fax Number: *For R & D use only*

Section 2 - Hazardous Ingredients/Identity Information

Emergency Overview OSHA Hazards

Flammable liquids (Category 2), H225 Acute toxicity, Oral (Category 4), H302 Acute toxicity, Inhalation (Category 4), H332 Acute toxicity, Dermal (Category 4), H312 Eye irritation (Category 2A), H319 GHS Label elements, including precautionary statements

Pictogram	$\forall \lor \lor \lor$	
Signal word	Danger	
Hazard statement(s)		
H225	Highly flammable liquid and vapour.	
H302 + H312 + H332	Harmful if swallowed, in contact with skin or if inhaled	
H319	Causes serious eye irritation.	
Precautionary statement(s		
P210	Keep away from heat/sparks/open flames/hot surfaces No smoking.	
P233	Keep container tightly closed.	
P240	Ground/bond container and receiving equipment.	
P241	Use explosion-proof electrical/ ventilating/ lighting/ equipment.	
P242	Use only non-sparking tools.Precautionary statement(s)	
P201	Obtain special instructions before use.	
P308 + P313	IF exposed or concerned: Get medical advice/ attention.	
P243	Take precautionary measures against static discharge.	
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.	
P264	Wash skin thoroughly after handling.	
P270	Do not eat, drink or smoke when using this product.	
P271	Use only outdoors or in a well-ventilated area.	
P280	Wear protective gloves/ eye protection/ face protection.	
P301 + P312 + P330 IF SW	ALLOWED: Call a POISON CENTER or doctor/ physician if you feel unwell. Rinse	
mouth.		
P303 + P361 + P353 IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin		
with water/ shower.		
P304 + P340 + P312 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for		
	breathing. Call a POISON	
	CENTER or doctor/ physician if you feel unwell.	
P305 + P351 + P338 IF IN	EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.	
P337 + P313	If eye irritation persists: Get medical advice/ attention.	

MSDS

Date prepared: July 25, 2014

Material Safety Data Sheet

Date revised: N/A (first edition)

P363	Wash contaminated clothing before reuse.
P370 + P378	In case of fire: Use dry sand, dry chemical or alcohol-resistant foam for extinction.
P403 + P235	Store in a well-ventilated place. Keep cool.
P501	Dispose of contents/ container to an approved waste disposal plant.

Section 3 - Hazardous Components/Identity Information

Hazardous Components	CAS #	Concentration
Acetonitrile	75-05-8	Balance
Melengestrol acetatesolution	2919-66-6	0.01%

Section 4 - First Aid Measures

General advice

Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area. If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician. In case of skin contact

Wash off with soap and plenty of water. Consult a physician.

In case of eye contact

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

If swallowed

Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

Section 5 - Fire-Fighting Measures

Extinguishing media

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special hazards arising from the substance or mixture

No data available

Advice for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

Further information

Use water spray to cool unopened containers.

Section 6 - Accidental Release Measures

Personal precautions

Use personal protective equipment. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Remove all sources of ignition. Evacuate personnel to safe areas. Beware of vapours accumulating to form explosive concentrations. Vapours can accumulate in low areas. Prevent further leakage or spillage if safe to do so. Do not let product

Environmental precautions enter drains. Methods for cleaning up

Contain spillage, and then collect with an electrically protected vacuum cleaner or by wet-brushing and place in container for disposal according to local regulations

Section 7 - Handling and Storage

Handling	Avoid contact with skin and eyes. Avoid inhalation of vapour or mist. Use	
explosion-proof		
	equipment.Keep away from sources of ignition - No smoking.Take measures to	
	prevent the build up of electrostatic charge.	
Storage	Keep container tightly closed in a dry and well-ventilated place. Containers	
	which are opened must be carefully resealed and kept upright to prevent	
	leakage.Store at -20°C.	

MSDS Date prepared: July 25, 2014

Material Safety Data Sheet

Date revised: N/A (first edition)

Section 8 - Exposure Controls/Personal Protection

Contains no substances with occupational exposure limit values

Personal protective equipment

Respiratory protection

Where risk assessment shows air-purifying respirators are appropriate use a full-face particle respirator type N100 (US) or type P3 (EN 143) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protectionFor prolonged or repeated contact use protective glovesEye protectionSafety glassesHygiene measuresGeneral industrial hygiene practice

Section 9 - Physical and Chemical Properties

Appearance:	Form	Liquid
	Color	Colorless
Safety data		
pH:		not available
Meltin	ng point:	not available
Boilin	g point:	not available
Flash	point:	not available
Ignitio	on temperature:	not available
Lowe	r explosion limit:	not available
Upper	explosion limit:	not available
Water	solubility:	not available
	3	

Section 10 - Stability and Reactivity

Storage stabilityStable under recommended storage conditionsMaterials to avoidHeat, flames and sparks. Extremes of temperature and direct sunlight.Hazardous decomposition productsHazardous decomposition products formed under fire
conditions. Carbon oxide, nitrogen oxides (NOx).

Section 11 - Toxicological Information

Acute toxicity:	no data available	
Other information on acute toxicity no data available		
Irritation and corrosion: not available		
Sensitization:	not available	
Respiratory or skin sensitization: not available		
Potential health effects		
Inhalation: May be harmful if inhaled. May cause respiratory tract irritation.		
Skin: May be harmful if absorbed through skin. May cause skin irritation.		
Eyes:May cause eye irritation.Ingestion:Toxic if swallowed.		

Section 12 - Ecological Information

Toxicity	not available	
Persistence and degradability	not available	
Bioaccumulative potential	not available	
Mobility in soil	not available	
Results of PBT and vPvB assessment not available		

MSDS

Material Safety Data Sheet

Date prepared: July 25, 2014

Date revised: N/A (first edition)

Section 13 - Disposal Considerations Product Burn in a ch

Burn in a chemical incinerator equipped with an afterburner and scrubber but exert extra care in igniting

as this material is highly flammable. Offer surplus and non-recyclable solutions to a licensed disposal

company. Contact a licensed professional waste disposal service to dispose of this material.

Contaminated packaging: Dispose of as unused product.

Section 14- Transport Information

\sim			
	DOT	IATA	
	Proper Shipping Name: Acetonitrile	Proper Shipping Name: Acetonitrile	
	UN#: 1648	UN#: 1648	
	Class: 3	Class: 3	
	Packing Group: II	Packing Group: II	
	Small quantities exempted	Small quantities exempted	

Section 15 - Regulatory Information

This product contains a chemical known to the State of California to cause birth defects or other reproductive harm.

Section 16 - Other Information

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. A Chemtek shall not be held liable for any damage resulting from handling or from contact with the above product.